


CLINICAL CENTRIFUGE



12 x 15ml Rotor



24 x 10ml Rotor



43mm Rubber Adapter



23mm Rubber Adapter



- :- Easy Programming
- :- Bright LCD Display
- :- 3 - Operations Mode
- :- Time & Speed Display Simultaneously
- :- RPM and RCF Both Can be displayed
- :- Noise Free Operation
- :- Sturdy & Compact Design
- :- Maintenance Free Brushless DC Motor
- :- Automatic Self diagnosis & Electric Lock for Better safety
- :- Stainless Steel Large Rotor for better Safety & Centrifugal force

Specifications	
Speed	300-4500 RPM
Max. RCF	Rotor 1 = 2558xg, Rotor 2 = 2200xg
Speed Accuracy	± 20 RPM
Timer	10s-99min
Motor	Brushless DC Motor
Operation Mode	Continuous / Short Time run / Timing
Safety Devices	Over Speed Detector, Door Interlock, Automatic Internal diagnosis
Rotor	12x15ml or 24x10ml
Noise Level	<65DB
Dimensions	450x340x230mm
Weight	10kg

WE **CARE** FOR
YOUR **RESEARCH** NEEDS

www.moxcare.in



Cat No.	Description
MCCLCG01	Clinical Centrifuge (Price Includes Any One Rotor)
MCCLCGRT1	Rotor 12x15ml Tube
MCCLCGRT2	Rotor 24x10ml Tube
RBA12	23mm/43mm Rubber Adapter 12 Pcs./Pack



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: info@moxcare.in W: www.moxcare.in