



- Ideal for tube sizes 1.5ml to 50ml
- Fixed Speed 15 RPM
- Optimized speed and tilt angle for thorough mixing in sample tubes
- There are no complicated controls, no setting or adjustment to be made just press the power switch and let it rock.

Gel Rocker

Specifications	
Motion	Rocking
Speed Range	15rpm
Maximum Loading Weight	2.5kg
Running Mode	Continuous
Rocking Angle	25°
Voltage / Power	AC 220V, 50/60Hz / 6W
Working plate size	252 x170mm
Dimensions	346 x 196 x 95mm
Net Weight	2.5kg

Cat No. MCGR1

WE **CARE** FOR
YOUR **RESEARCH** NEEDS

www.moxcare.in



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: info@moxcare.in W: www.moxcare.in