

## MAGNETIC STIRRER

Compact and Sturdy design



## STIRRMAG mini

- Compact and Sturdy design
- 1Lt. Stirring Capacity
- Burshless DC Motor
- High Temp. Resistant Plastic Body

Specifications	STIRRMAG mini
Operation Mode	Stirrer
Speed Range	100-1500RPM
Max Stirring Capacity (H2O)	1000ml
Dimensions of Plate	120mm
Working Plate Material	Plastic
Display	n/a
Temperature Range	n/a
Temperature Display Accuracy	n/a
External Temperature Sensor	n/a
External Temp.Sensor Accuracy	n/a
Timer	n/a
Hot Surface Warning	n/a
Heating Output	n/a
Motor Type	Brushless DC Motor
Permissible Temp./Humidity	5-40°C, 80%
Protection class	Ip21
Voltage, Frequency	200-240V, 50/60Hz
Power	10W
Dimensions	138x60mm
Weight	0.7Kg

Catalog No.	Description
ISTIRRMAG mini	Mini Magnetic Stirrer
MCSTRNG AP	Stirring bars Assorted Pack 15x8(2Pcs.),20x8 (2Pcs.),25x8(2Pcs.), 40x8 (1Pcs.), 50x8 (1Pcs.)
MCSTRNGX	Stirring Bar Cross Shape (2pcs./Pkg)
MCSTRNGQ	Stirring Bar Quincunx Shape (2pcs./Pkg)

WE **CARE** FOR  
YOUR **RESEARCH** NEEDS

[www.moxcare.in](http://www.moxcare.in)



## MOXCARE PRODUCTS INC.

Unit 1  
115, HSIIDC Industrial Area  
Ambala Cantt - 133006 Haryana-India

Unit 2  
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center  
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829  
E: [info@moxcare.in](mailto:info@moxcare.in) W: [www.moxcare.in](http://www.moxcare.in)