

Minimax10000 MINI CENTRIFUGE



- Variable Speed 3000-10000RPM
- Low Noise ≤ 50 dB
- Enclosed Design Reduces risk of liquid splashing into equipment
- Click-On Design allows rotor changing very Easy
- Equipped with 2 interchangeable rotors 8 x 0.5ml/1.5ml/2ml centrifuge tube or 2 x 0.2ml Strip of 8

| Specifications | Minimax 10000 |
|-----------------|---|
| Speed Range | 3000-10000 RPM (1000rpm/grade) |
| Maximum RCF | 5460xg |
| Running Mode | Timing/Continuous |
| Timer | 1s ~ 99min |
| Rotors Capacity | 8x0.5ml/1.5ml/2ml centrifuge tube or 2 x 0.2ml Strip of 8 |
| Noise Range | <50dB |
| Voltage | AC 220V, 50/60Hz |
| Dimensions | 195(L) x 170(W) x 135(H) |
| Net Weight | 0.8kg |

Cat No. MINIMAX10000

WE **CARE** FOR
YOUR **RESEARCH** NEEDS

www.moxcare.in



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: info@moxcare.in W: www.moxcare.in