

High Speed Mini Centrifuge

Mini 15000



- Brushless DC Motor
- Large and Bright LED Display
- Powerful centrifugal force 15100 x g and Max.Speed 15000 RPM
- PULSE to offer Quick Spins Explosion Proof Inner Steel Build Design
- Heavy Base Eliminates Vibration & Movement when operating at High Speeds
- Fast Acceleration & Braking Digital Wheel to adjust Speed and Time
- Motor Driven Lid-Lock System for Safety & reliability
- Programmable - Can store upto 10 Programs
- Automatic Lid Open upon running completion
- *Timer Starts when machine reach the set speed

Specifications	
Speed	Variable 500-15000RPM,
Max.RCF [x g]	15100g
Increment	100RPM / 100 RCF
Speed accuracy [rpm]	±20
Rotor capacity	0.2ml/0.5ml/1.5ml/2ml x 12
Timer	30sec-99min
Mode	Continuous/Timer
Motor	Brushless DC motor
Noise level [db]	≤56
Acceleration/Braking time [Sec]	30sec / 10sec
Safety devices	Over-Speed detection, Over-Temperature detection Door interlock & Auto Internal Diagnosis
Additional features	Speed RCF Switch, Short-Time Run Function, Sound-Alert function
Power	Single-phase, AC110V-240V-,50/Hz/60Hz,3A
Dimensions [mm]	256 x 225 x 155
Weight	6kg

Cat No.	Description
MINI15000+	High Speed Mini Centrifuge
M15RTK	Rotor Adapter Kit for 0.2ml



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: info@moxcare.in W: www.moxcare.in