


Minimax5400 MINI CENTRIFUGE



- Fixed Speed 5400 RPM
- Low Noise \leq 50dB
- Enclosed Design Reduces risk of liquid splashing into equipment
- Click-On Design allows rotor changing very Easy
- Equipped with 2 interchangeable rotors 1.5/2ml x 6 and PCR 8x2

Specifications	
Speed	5400 RPM
Max RCF [x g]	2000g
Motor	DC motor
Noise level	≤50
Mode	Continuous
Rotor	0.2/0.5/1.5/2.0ml x 6 & 0.2ml x 2 PCR 8 strips
Input	100-240VAC,50/60 HZ
Output	24VDC,2A
Dimension	6.9" x 5.5"
Weight	1.85kg

Cat No. MINIMAX5400

WE **CARE** FOR
YOUR **RESEARCH** NEEDS

www.moxcare.in



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: info@moxcare.in W: www.moxcare.in