



Specifications	
Speed Range	50-300 RPM
Shaking Orbit	20mm
Temperature Range	RT+5°C-60°C
Display Accuracy	±0.1°C
Timer	1min-999min
Display	Super Bright LED
Controller	Programmable Microprocessor
Power	400W
Voltage	220V, 50/60Hz
Dimensions	360x435x320mm
Weight	18.7kg

Cat No. Virtue ISC300

- Compact size
- Excellent Temperature uniformity
- Built-in Cover Switch
- Perfect for low capacity application
- Brushless DC motor ensures smooth operation
- Beeps when Program Completes
- Temperature Protection



MOXCARE PRODUCTS INC.

Unit 1
115, HSIIDC Industrial Area
Ambala Cantt - 133006 Haryana-India

Unit 2
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center
Saha, Ambala - 133104 Haryana-India


T: +91 171 2699017, 829 5958 829
E: info@moxcare.in W: www.moxcare.in



MCOS-UPAB



MCOS-C10250/
MCOS-C25300



MCOS-FP