



# WATER BATH

- Compact unique Design
- PID automatic temperature control
- Built in water level sensor
- Separate Display for temperature and time
- Easy user interface

Specifications	Water Bath 2L	Water Bath 5L
Display	LED Display	
Temperature Range	RT.+5 - 100°C	
Temperature Accuracy	0.3°C	
Temperature Stability at 37°C	0.5°C	
Timer	Yes	
Both Volume	2L x 1	5L x 1
Inner Dimension	135mm x 150mm x 143mm x 1	325mm x 150mm x 176mm x 1
Parts Touch Water	Stainless Steel	
Voltage / Power	AC 220V 50/60Hz /300W	AC 220V, 50/60Hz /1000W
Dimension / Net Weight	198mm x 290mm 240mm / 2.3kg.	360mm x 286mm x 245mm / 5.1kg

Cat No. MCWB02,MCWB05

WE **CARE** FOR  
YOUR **RESEARCH** NEEDS

[www.moxcare.in](http://www.moxcare.in)



## MOXCARE PRODUCTS INC.

Unit 1  
115, HSIIDC Industrial Area  
Ambala Cantt - 133006 Haryana-India

Unit 2  
308-E, Sec-2, Phase-1, HSIIDC Industrial Growth Center  
Saha, Ambala - 133104 Haryana-India

T: +91 171 2699017, 829 5958 829

E: [info@moxcare.in](mailto:info@moxcare.in) W: [www.moxcare.in](http://www.moxcare.in)